







LABORATORY MINERAL JIG



Single Source

Incorporated testing, design engineering, fabrication, and services to confirm tailor-made total systems

Star Trace offers a range of products which are widely used across various laboratories and also hold diverse industrial applications. These wide range of products are manufactured according to industrial quality and standards. Backed by our excellent engineering and quality control team, we are able to satisfy the client requirement for various engineering products all over the world. With more than 20 years of experience in this field, we have been supplying technically accurate and superior products to our clients and have catered to the requirement of variety of sectors and clients.

MINERAL JIG

Principle: Jigging is the process of separating the particles of different specific gravity, size and shape by introducing them on screen though which a fluid is made to pulsate alternately. It is consists of a rectangular open top container. The bottom of the container is either a screen or a perforated surface which supports the bed of the solids. Pulsion and suction strocks of water through the screen are given by means of a reciprocating piston.

During the pulsion stroke, solid bed is opened and the particles are lifted. The particles settle into more homogeneous layer during the suction stroke. By repeated pulsion and suction stroke, a complete stratification is achieved with the heavy particles at the bottom and lighter particles above them.

The following three effects contributes to the stratification in Jig:-

- 1. Hindered settling classification.
- 2. Differential acceleration at the beginning of fall.
- 3. Consolidation trickling at the end of the fall.

Model	Area m2	Stroke mm	Times of Stroke r/min	Feed Size	Capacity t/h	Feedwater Flow Control m3/h	Power kw	Overall Dimension mm	Weight t
JT2-2	2.28	12-17	80-180	<10	4-8	1-3	3	3220/1550/2050	1.6
JT3-1	3.3		80-120	<10	6-10	1-3	5.5	3180/1850/2600	3.26
JT4-2	4.33	6.5-25	80-120	<30	10-15	1-2	7.5	3600/1850/2600	4.6
JT5-2	4.86	15-25	80-120	<10	10-20	1-2	7.5	3600/2000/26000	4.8
LTA1010/2	4	0-26	250-350	1-6	2-20	1-2	2.2	2800/1800/2300	1.3
300/450 Diaphragm	2.7	0-26	320	12	3-6	2-3	1.1	1230/1270/1750	0.74
370/360 Diaphragm	2.75	5 25	2-5	6	1-3	2-3	1.1	1000/6000/1600	0.24

Single Source

Corporate Headquarters

Old No: 7, New No: 5, Jeevanandham Street, Redhills, Chennai-600 052, India. Phone: (0091-044) 26418456, 26419439, 32975456. TeleFax: (0091-044) 26419439 Mobile: (0091) 9380521576, 9380521567, 9841015456. Email: sales@startraceltd.com

City Office

Old No: 526, New No: 41, 1" Floor, 3rd Avenue Road, N.S.K Nagar, Arumbakkam, Chennal-600 106, India, TeleFax: (0091-044) 43500032

Works

lo:11, Moorthasa Ussian Sahib Street, Redhills, Chennai-600 052, India.

Canadian Representative

Mr. Mahes parteepan Phone: (001) 2988003549 Mobile: (001) 647 783 9627 Email: mahesh@startraceltd.co

Myanmar Representative Mr. KYAW MIN TUN

Mobile: +95-973166695 Email: kyaw@startrace.in

Indonesia Office

BLOK D6-N0.5, MGK - Mega Kemayoran, Office Tower A Ll.1 Jl. Angkasa, Kota Baru Bandar Kemayoran, Jakarta Pusat 10610. Phone: +62 81311594498. Fmail'sales@startrace co.id

Dubai Representative

Mobile: 00971557526778 Fmail: mahesh@startraceltd.com